



Use your own control switches

Big Jack - Remote control

ART NR: 402500

Big Jack makes it easy for you to e.g. watch your favourite channel on TV, listen to the radio, switch on the light, open the door to visitors, or answer the telephone.

Big Jack is a robust, programmable remote control with several optional, external control switches that allow you to easily manage various media devices and functions in your environment. It has been specially designed for people who suffer from cognitive or motor skill disabilities, and who therefore require large and clear control buttons.

Control buttons are available in several different colours, and you are able to label them with pictures, symbols or text. You are able to clearly see which device or function you are controlling by using that exact button.



ABILIA

MANAGE EVERYTHING USING EXTERNAL CONTROL BUTTONS

Big Jack is a robust remote control that you can use to manage your media devices and various functions in your home.

The large, clear, outer control buttons make it easy for you to use your TV, stereo, radio, loudspeaker, CD-player, or to control your lights, doors, windows and telephone.

A REMOTE CONTROL THAT CAN HANDLE ALL YOUR NEEDS

Big Jack is programmable and can “learn” functions via IR signals from your other remote controls, such as remote controls for the TV, loudspeaker system and stereo. That means that you only need to use one remote control to manage all your devices and functions.

AUTOMATICALLY CONTROL A SERIES OF EVENTS

You can create sequences that allow a single button to replace an entire series of button pushes, and which automatically start a series of events. This would, for example, make it easier to find a specific teletext channel, your favourite TV show on SVT-play, or a film channel.

An additional application could be that the sound on the TV switches off when you answer your door's intercom.

COMES WITH

User manual.
Batteries and control buttons must be ordered separately.



EASY TO USE

You can connect up to six different control buttons on your Big Jack. They are available in several different colours, and you can label them with pictures, symbols or text to simply show exactly what each button does.

It could not get simpler than that.



TECHNICAL INFORMATION		PRODUCT	ART NR
Size	140 x 100 x 38 mm	Big Jack IR-6SP	402500
Weight	280 g		
Coverage	20 - 30 m	ACCESSORIES	
Input sockets	6 x 3,5 mm	Pikobutton 30 mm, Red	451000
		Pikobutton 50 mm, Red	451500
		Extension cable 3,5 mm	4538

Ver. 2023.02.03